


DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

SECTION 1 - PRODUCT IDENTIFICATION AND USE/IDENTIFICATION ET UTILISATION DU PRODUIT				
PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT		AH 26 Powder Silver Free		
PRODUCT USE USAGE DU PRODUIT		Dental		
MANUFACTURER'S NAME/NOM DU FABRICANT		SUPPLIER'S NAME/NOM DU FOURNISSEUR DENTSPLY Canada Ltd.		
STREET ADDRESS/NUMERO CIVIQUE De-Tray-Str.1		STREET ADDRESS/NUMERO CIVIQUE 161 Vinyl Court		
CITY/VILLE Konstanz	PROVINCE Germany	CITY/VILLE Woodbridge	PROVINCE Ontario	
POSTAL CODE POSTAL D-78467	EMERGENCY TELEPHONE NO. D'URGENCE 1-49-7531-583-0	POSTAL CODE POSTAL L4L 4A3	EMERGENCY TELEPHONE NO. D'URGENCE (905) 851-6060 or 1-800-263-1437 FAX# (905) 851-9809	
SECTION 2 - HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX				
HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX	%	CAS NUMBER/ NUMERO CAS	TLV	OSHA PEL
Methenamine - Hexamethylenetetraamine	20	100-97-0		
SECTION 3 - PHYSICAL DATA/CHARACTERISTIQUES PHYSIQUES				
PHYSICAL STATE/ETAT PHYSIQUE Powder	ODOUR AND APPEARANCE/ODEUR ET APPARENCE: Yellow powder.		Chemical Family: Amine.	
VAPOUR PRESSURE (mm Hg)/ TENSION DE VAPEUR (mm Hg) N.A.	VAPOUR DENSITY (AIR-1)/ DENSITE DE LA VAPEUR (AIR-1) N.A.	EVAPORATION RATE/ TAUX D'EVAPORATION N.A.	BOILING POINT (°C)/ POINT D'EBULLITION (°C) N.A.	FREEZING POINT (°C) POINT DE CONGELATION (°C) N.A.
pH: N.A.	SPECIFIC GRAVITY/DENSITE: N.A.	SOLUBILITY IN WATER: Partly Soluble.		
SECTION 4 - FIRE AND EXPLOSION DATA/RISQUES D'INCENDIE OU D'EXPLOSION				
FLAMMABILITY/INFLAMMABILITE				
YES/OUI <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF YES, UNDER WHICH CONDITIONS? At elevated temperatures highly flammable, fine powder can present dust explosion hazard. SI OUI, DANS QUELLES CONDITIONS?				
MEANS OF EXTINCTION/MOYENS D'EXTINCTION: Water spray, Foam, CO ₂				
FLASHPOINT (°C) AND METHOD/ POINT D'ECLAIR (°C) ET METHODE DE DETERMINATION 250°C		UPPER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MAXIMAL (% PAR VOLUME) N.A.	LOWER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MINIMAL (% PAR VOLUME) N.A.	
AUTO-IGNITION TEMPERATURE (°C)/ TEMPERATURE D'AUTO-INFLAMMATION (°C): N.A.		HAZARDOUS COMBUSTION PRODUCTS/ PRODUITS DANGEREUX DE COMBUSTION: Poisonous gases may be released during combustion which can form explosive gas-air mixture.		
EXPLOSION DATA/ DONNES SUR LES EXPLOSIONS		SENSITIVITY TO STATIC DISCHARGE SENSIBILITE A UNE DECHARGE STATIQUE N.A.		
SECTION 5 - REACTIVITY DATA/DONNES SUR LA REACTIVITE				
CHEMICAL STABILITY/STABILITE CHIMIQUE				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF NO, UNDER WHICH CONDITIONS? SI NON, DANS QUELLES CONDITIONS?				
INCOMPATIBILITY WITH OTHER SUBSTANCES/ INCOMPATIBILITE AVEC D'AUTRES PRODUITS				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF SO, WHICH ONES? Oxidizing chemicals. SI OUI, LESQUELS?				
REACTIVITY, AND UNDER WHAT CONDITIONS?/CONDITIONS DE REACTIVITE Elevated temperatures.				
HAZARDOUS DECOMPOSITION PRODUCTS/PRODUITS DE DECOMPOSITION DANGEREUX Formaldehyde, Nitrogen oxides, Ammonia				

DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT				AH 26 Powder Silver Free			
SECTION 6 - TOXICOLOGICAL PROPERTIES/PROPRIÉTÉS TOXICOLOGIQUES							
ROUTE OF ENTRY/VOIES D'ABSORPTION							
SKIN CONTACT/ CONTACT DERMIQUE		SKIN ABSORPTION/ ABSORPTION DERMIQUE		EYE CONTACT/ CONTACT OCULAIRE		INHALATION	INGESTION
X		—		X		—	X
EFFECTS OF ACUTE EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION AIGUE AU PRODUIT							
N.A.							
EFFECTS OF CHRONIC EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION CHRONIQUE AU PRODUIT							
EXPOSURE LIMITS/ SEUILS LIMITES D'EXPOSITION							
N.A.		IRRITANCY OF PRODUCT/ PROPRIÉTÉ IRRITANTE		SENSITIZATION TO PRODUCT/ SENSIBILISATION AU PRODUIT		CARCINOGENICITY/ CANCÉROGÉNÉCITÉ	
		Yes		Yes		N.A.	
TERATOGENICITY/ TERATOGENICITÉ		REPRODUCTIVE TOXICITY/ EFFETS SUR LA REPRODUCTION		MUTAGENICITY/ MUTAGÉNÉCITÉ		SYNERGISTIC PRODUCTS/ PRODUITS SYNERGIQUES	
N.A.		N.A.		N.A.		N.A.	
SECTION 7 - PREVENTIVE MEASURES/MESURES PREVENTIVES							
PERSONAL PROTECTIVE EQUIPMENT/MATÉRIEL PERSONNEL DE PROTECTION							
Use protective gloves and safety glasses.							
GLOVES (SPECIFY)/GANTS (PRÉCISER):		RESPIRATOR (SPECIFY)/APPAREIL RESPIRATOIRE (PRÉCISER)		EYE (SPECIFY)/APPAREIL OCULAIRE (PRÉCISER):			
Impervious Gloves.		N.A.		Protective Goggles.			
FOOTWEAR (SPECIFY)/CHAUSSURES (PRÉCISER)		CLOTHING (SPECIFY)/VÊTEMENTS (PRÉCISER)		OTHER (SPECIFY)/AUTRES (PRÉCISER)			
N.A.		N.A.					
ENGINEERING CONTROLS (SPECIFY, EG, VENTILATION, ENCLOSED PROCESS)/INSTALLATION TECHNIQUES							
Ensure good ventilation and local exhaust. Prevent formation of dust.							
LEAK AND SPILL PROCEDURE/PROCÉDURES EN CAS DE FUITES ET DE DEVERSEMENTS							
Eliminate ignition sources. Avoid breathing dust.							
WASTE DISPOSAL/ÉLIMINATION DES RÉSIDUS:							
Disposal must be made according to the local regulations. Do not dispose in municipal waste. Do not allow to reach sewage system.							
HANDLING PROCEDURES AND EQUIPMENT/MÉTHODES ET ÉQUIPEMENT POUR LA MANUTENTION:							
Store at ambient temperatures out of direct sunlight. Use adequate ventilation.							
STORAGE REQUIREMENTS/EXIGENCES EN MATIÈRE D'ENTREPOSAGE:							
Avoid formation of dust. Do not allow to reach sewer system or water course. Ensure good ventilation.							
SPECIAL SHIPPING INFORMATION/RENSEIGNEMENTS SPÉCIAUX EN MATIÈRE D'EXPÉDITION							
Flammable Solid UN 1328, Class 4.1, Packing Group III.							
SECTION 8 - FIRST AID MEASURES/PREMIER SOINS							
SPECIFIC MEASURES/PREMIERS SOINS PARTICULIERS À ADMINISTRER							
Eyes: Flush immediately with water for minimum of 15 minutes and seek medical attention.							
Skin: (Irritation). Wash with soap and water, remove contaminated clothing.							
Inhalation: Move subject to fresh air.							
Ingestion: Induce vomiting and consult physician immediately.							
SECTION 9 - PREPARATION DATE OF MSDS/RENSEIGNEMENTS SUR LA PRÉPARATION DE LA FS							
PREPARED BY (GROUP, DEPARTMENT, ETC.) PRÉPARÉE PAR (GROUPE, SERVICE, PERSONNE, ETC.)		PHONE NUMBER/NO. DE TÉLÉPHONE		DATE			
 Cadentrix		(416) 994-9356		Rev. 2 1 June 2005			

The statements contained herein are offered for informational purposes only and are based upon technical data that Dentsply Ltd. believes to be accurate. It is intended for use only by persons having the necessary technical skill and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, and FITNESS OR OTHERWISE.