


DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

SECTION 1 - PRODUCT IDENTIFICATION AND USE/IDENTIFICATION ET UTILISATION DU PRODUIT				
PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT		AH Plus™ Paste A		
		PRODUCT IDENTIFICATION NUMBER (PIN)/ NUMERO D'IDENTIFICATION DU PRODUIT (NIP) 606 20110		
PRODUCT USE USAGE DU PRODUIT Dental				
MANUFACTURER'S NAME/NOM DU FABRICANT DENTSPLY DeTray GmbH			SUPPLIER'S NAME/NOM DU FOURNISSEUR DENTSPLY Canada Ltd.	
STREET ADDRESS/NUMERO CIVIQUE De-Tray-Str.1			STREET ADDRESS/NUMERO CIVIQUE 161 Vinyl Court	
CITY/VILLE Konstanz	PROVINCE Germany	CITY/VILLE Woodbridge	PROVINCE Ontario	
POSTAL CODE POSTAL D-78467	EMERGENCY TELEPHONE NO. D'URGENCE 1-49-7531-583-0	POSTAL CODE POSTAL L4L 4A3	EMERGENCY TELEPHONE NO. D'URGENCE (905) 851-6060 or 1-800-263-1437 FAX# (905) 851-9809	
SECTION 2 - HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX				
HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX	%	CAS NUMBER/ NUMERO CAS	TLV	OSHA/PEL
Epoxide resin	30	N.A.	LD₅₀ > 5000 mg/Kg	
SECTION 3 - PHYSICAL DATA/CHARACTERISTIQUES PHYSIQUES				
PHYSICAL STATE/ETAT PHYSIQUE	ODOUR AND APPEARANCE/ODEUR ET APPARENCE:		Chemical Family:	
Paste	Odourless Yellow paste.		Epoxy	
VAPOUR PRESSURE (mm Hg)/ TENSION DE VAPEUR (mm Hg)	VAPOUR DENSITY (AIR-1)/ DENSITE DE LA VAPEUR (AIR-1)	EVAPORATION RATE/ TAUX D'EVAPORATION	BOILING POINT (°C)/ POINT D'EBULLITION (°C)	FREEZING POINT (°C) POINT DE CONGELATION (°C)
N.A.	N.A.	N.A.	N.A.	
pH:	SPECIFIC GRAVITY/DENSITE:	SOLUBILITY IN WATER:		
N.A.	N.D.	Insoluble.		
SECTION 4 - FIRE AND EXPLOSION DATA/RISQUES D'INCENDIE OU D'EXPLOSION				
FLAMMABILITY/INFLAMMABILITE				
YES/OUI <input type="checkbox"/> NO/NON <input checked="" type="checkbox"/> IF YES, UNDER WHICH CONDITIONS? SI OUI, DANS QUELLES CONDITIONS?				
MEANS OF EXTINCTION/MOYENS D'EXTINCTION: Water spray, Foam, CO₂				
FLASHPOINT (°C) AND METHOD/ POINT D'ECLAIR (°C) ET METHODE DE DETERMINATION	UPPER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MAXIMAL (% PAR VOLUME)		LOWER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MINIMAL (% PAR VOLUME)	
N.A.	N.A.		N.A.	
AUTO-IGNITION TEMPERATURE (°C)/ TEMPERATURE D'AUTO-INFLAMMATION (°C):	HAZARDOUS COMBUSTION PRODUCTS/ PRODUITS DANGEREUX DE COMBUSTION: Irritating gases may be released during combustion.			
N.A.				
SENSITIVITY TO IMPACT/SENSIBILITE A L'IMPACT MECANIQUE			SENSITIVITY TO STATIC DISCHARGE SENSIBILITE A UNE DECHARGE STATIQUE	
EXPLOSION DATA/ DONNES SUR LES EXPLOSIONS N.A.			N.A.	
SECTION 5 - REACTIVITY DATA/DONNES SUR LA REACTIVITE				
CHEMICAL STABILITY/STABILITE CHIMIQUE				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF NO, UNDER WHICH CONDITIONS? SI NON, DANS QUELLES CONDITIONS?				
INCOMPATIBILITY WITH OTHER SUBSTANCES/ INCOMPATIBILITE AVEC D'AUTRES PRODUITS				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF SO, WHICH ONES? Alkali metals, oxidizing chemicals, Halogens. SI OUI, LESQUELS?				
REACTIVITY, AND UNDER WHAT CONDITIONS?/CONDITIONS DE REACTIVITE				
Avoid high temperatures and humidity.				
HAZARDOUS DECOMPOSITION PRODUCTS/PRODUITS DE DECOMPOSITION DANGEREUX				
Irritating gases.				

DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT			
AH Plus™ Paste A			
SECTION 6 - TOXICOLOGICAL PROPERTIES/PROPRIÉTÉS TOXICOLOGIQUES			
ROUTE OF ENTRY/VOIES D'ABSORPTION			
SKIN CONTACT/ CONTACT DERMIQUE	X	SKIN ABSORPTION/ ABSORPTION DERMIQUE	—
		EYE CONTACT/ CONTACT OCULAIRE	X
		INHALATION	—
		INGESTION	—
EFFECTS OF ACUTE EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION AIGUE AU PRODUIT Material is non-toxic. Contact dermatitis (Skin sensitization) can develop in sensitive individuals. Avoid direct contact with skin.			
EFFECTS OF CHRONIC EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION CHRONIQUE AU PRODUIT Material is non-toxic. Contact dermatitis (Skin sensitization) can develop in sensitive individuals.			
EXPOSURE LIMITS/ SEUILS LIMITES D'EXPOSITION N.A.	IRRITANCY OF PRODUCT/ PROPRIÉTÉ IRRITANTE Yes	SENSITIZATION TO PRODUCT/ SENSIBILISATION AU PRODUIT Yes	CARCINOGENICITY/ CANCÉROGENICITÉ N.A.
TERATOGENICITY/ TERATOGENICITÉ N.A.	REPRODUCTIVE TOXICITY/ EFFETS SUR LA REPRODUCTION N.A.	MUTAGENICITY/ MUTAGÉNICITÉ N.A.	SYNERGISTIC PRODUCTS/ PRODUITS SYNERGIQUES N.A.
SECTION 7 - PREVENTIVE MEASURES/MESURES PREVENTIVES			
PERSONAL PROTECTIVE EQUIPMENT/MATÉRIEL PERSONNEL DE PROTECTION			
GLOVES (SPECIFY)/GANTS (PRÉCISER): Users with special skin problems, cuts or abrasions should wear protective gloves (Latex or PVC).	RESPIRATOR (SPECIFY)/APPAREIL RESPIRATOIRE (PRÉCISER) When working with large amount wear a mask.	EYE (SPECIFY)/APPAREIL OCULAIRE (PRÉCISER): Protective Goggles.	
FOOTWEAR (SPECIFY)/CHAUSSURES (PRÉCISER) N.A.	CLOTHING (SPECIFY)/VÊTEMENTS (PRÉCISER) Use protective clothing,	OTHER (SPECIFY)/AUTRES (PRÉCISER)	
ENGINEERING CONTROLS (SPECIFY, EG, VENTILATION, ENCLOSED PROCESS)/INSTALLATION TECHNIQUES			
LEAK AND SPILL PROCEDURE/PROCÉDURES EN CAS DE FUITES ET DE DEVERSEMENTS Absorb using universal absorbing material.			
WASTE DISPOSAL/ÉLIMINATION DES RÉSIDUS: Do not empty into drains, water or ground. Dispose according to local regulations.			
HANDLING PROCEDURES AND EQUIPMENT/MÉTHODES ET ÉQUIPEMENT POUR LA MANUTENTION: Refrigerated storage prolongs shelf life.			
STORAGE REQUIREMENTS/EXIGENCES EN MATIÈRE D'ENTREPOSAGE: Store away from food and beverages. Do not store near high temperature or ignition sources.			
SPECIAL SHIPPING INFORMATION/RENSEIGNEMENTS SPÉCIAUX EN MATIÈRE D'EXPÉDITION Not Regulated.			
SECTION 8 - FIRST AID MEASURES/PREMIER SOINS			
SPECIFIC MEASURES/PREMIERS SOINS PARTICULIERS À ADMINISTRER Eyes: Flush immediately with water for minimum of 15 minutes and seek medical attention. Skin: (Irritation). Wash with soap and water, remove contaminated clothing.			
SECTION 9 - PREPARATION DATE OF MSDS/RENSEIGNEMENTS SUR LA PRÉPARATION DE LA FS			
PREPARED BY (GROUP, DEPARTMENT, ETC.) PRÉPARÉ PAR (GROUPE, SERVICE, PERSONNE, ETC.)  Cadentrix	PHONE NUMBER/NO. DE TÉLÉPHONE (416) 994-9356	DATE Rev. 2 1 June 2005	

The statements contained herein are offered for informational purposes only and are based upon technical data that Dentsply Ltd. believes to be accurate. It is intended for use only by persons having the necessary technical skill and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, and FITNESS OR OTHERWISE.