


## DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

SECTION 1 - PRODUCT IDENTIFICATION AND USE/IDENTIFICATION ET UTILISATION DU PRODUIT				
PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT		AH Plus™ Paste B		
PRODUCT USE USAGE DU PRODUIT		Dental		
MANUFACTURER'S NAME/NOM DU FABRICANT		SUPPLIER'S NAME/NOM DU FOURNISSEUR		
DENTSPLY DeTray GmbH		DENTSPLY Canada Ltd.		
STREET ADDRESS/NUMERO CIVIQUE		STREET ADDRESS/NUMERO CIVIQUE		
De-Tray-Str.1		161 Vinyl Court		
CITY/VILLE	PROVINCE	CITY/VILLE	PROVINCE	
Konstanz	Germany	Woodbridge	Ontario	
POSTAL CODE POSTAL	EMERGENCY TELEPHONE NO. D'URGENCE	POSTAL CODE POSTAL	EMERGENCY TELEPHONE NO. D'URGENCE	
D-78467	1-49-7531-583-0	L4L 4A3	(905) 851-6060 or 1-800-263-1437 FAX# (905) 851-9809	
SECTION 2 - HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX				
HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX	%	CAS NUMBER/ NUMERO CAS		
Adamantine Amine	5	768-94-5	LD <sub>50</sub> 900 mg/Kg	
N,N-Dibenzyl-5 Oxanonanediamine-1,9	10	113506-22-2	LD <sub>50</sub> 4581 G/Kg	
TCD - Diamine	<1	68889-71-4	LD <sub>50</sub> 502 mg/Kg	
SECTION 3 - PHYSICAL DATA/CHARACTERISTIQUES PHYSIQUES				
PHYSICAL STATE/ETAT PHYSIQUE	ODOUR AND APPEARANCE/ODEUR ET APPARENCE:		Chemical Family:	
Paste	White paste with characteristic odour.		Amine.	
VAPOUR PRESSURE (mm Hg)/ TENSION DE VAPEUR (mm Hg)	VAPOUR DENSITY (AIR-1)/ DENSITE DE LA VAPEUR (AIR-1)	EVAPORATION RATE/ TAUX D'EVAPORATION	BOILING POINT (°C)/ POINT D'EBULLITION (°C)	FREEZING POINT (°C)/ POINT DE CONGELATION (°C)
N.A.	N.A.	N.A.	N.A.	N.A.
pH:	SPECIFIC GRAVITY/DENSITE:	SOLUBILITY IN WATER:		
N.A.	N.D.	Negligible.		
SECTION 4 - FIRE AND EXPLOSION DATA/RISQUES D'INCENDIE OU D'EXPLOSION				
FLAMMABILITY/INFLAMMABILITE				
YES/OUI <input type="checkbox"/> NO/NON <input checked="" type="checkbox"/>		IF YES, UNDER WHICH CONDITIONS? SI OUI, DANS QUELLES CONDITIONS?		
MEANS OF EXTINCTION/MOYENS D'EXTINCTION: Water spray, Foam, CO <sub>2</sub>				
FLASHPOINT (°C) AND METHOD/ POINT D'ECLAIR (°C) ET METHODE DE DETERMINATION	UPPER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MAXIMAL (% PAR VOLUME)		LOWER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MINIMAL (% PAR VOLUME)	
N.A.	N.A.		N.A.	
AUTO-IGNITION TEMPERATURE (°C)/ TEMPERATURE D'AUTO-INFLAMMATION (°C):	HAZARDOUS COMBUSTION PRODUCTS/ PRODUITS DANGEREUX DE COMBUSTION: Irritating gases may be released during combustion.			
N.A.				
EXPLOSION DATA DONNES SUR LES EXPLOSIONS		SENSITIVITY TO IMPACT/SENSIBILITE A L'IMPACT MECANIQUE		SENSITIVITY TO STATIC DISCHARGE SENSIBILITE A UNE DECHARGE STATIQUE
N.A.		N.A.		N.A.
SECTION 5 - REACTIVITY DATA/DONNES SUR LA REACTIVITE				
CHEMICAL STABILITY/STABILITE CHIMIQUE				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/>		IF NO, UNDER WHICH CONDITIONS? SI NON, DANS QUELLES CONDITIONS?		
INCOMPATIBILITY WITH OTHER SUBSTANCES/ INCOMPATIBILITE AVEC D'AUTRES PRODUITS				
YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/>		IF SO, WHICH ONES? SI OUI, LESQUELS? Alkali metals, oxidizing chemicals, Halogens.		
REACTIVITY, AND UNDER WHAT CONDITIONS?/CONDITIONS DE REACTIVITE				
None.				
HAZARDOUS DECOMPOSITION PRODUCTS/PRODUITS DE DECOMPOSITION DANGEREUX				
Irritating and toxic gases.				

## DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT			
<b>AH Plus™ Paste B</b>			
<b>SECTION 6 - TOXICOLOGICAL PROPERTIES/PROPRIÉTÉS TOXICOLOGIQUES</b>			
ROUTE OF ENTRY/VOIES D'ABSORPTION			
SKIN CONTACT/ CONTACT DERMIQUE	<b>X</b>	SKIN ABSORPTION/ ABSORPTION DERMIQUE	—
		EYE CONTACT/ CONTACT OCULAIRE	<b>X</b>
		INHALATION	—
		INGESTION	—
EFFECTS OF ACUTE EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION AIGUE AU PRODUIT <b>Material is non-toxic. Contact dermatitis (Skin sensitization) can develop in sensitive individuals.</b>			
EFFECTS OF CHRONIC EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION CHRONIQUE AU PRODUIT <b>Material is non-toxic. Contact dermatitis (Skin sensitization) can develop in sensitive individuals.</b>			
EXPOSURE LIMITS/ SEUILS LIMITES D'EXPOSITION	<b>N.A.</b>	IRRITANCY OF PRODUCT/ PROPRIÉTÉ IRRITANTE	<b>Yes</b>
		SENSITIZATION TO PRODUCT/ SENSIBILISATION AU PRODUIT	<b>Yes</b>
		CARCINOGENICITY/ CANCÉROGENICITÉ	<b>N.A.</b>
TERATOGENICITY/ TERATOGENICITÉ	<b>N.A.</b>	REPRODUCTIVE TOXICITY/ EFFETS SUR LA REPRODUCTION	<b>N.A.</b>
		MUTAGENICITY/ MUTAGÉNICITÉ	<b>N.A.</b>
		SYNERGISTIC PRODUCTS/ PRODUITS SYNERGIQUES	<b>N.A.</b>
<b>SECTION 7 - PREVENTIVE MEASURES/MESURES PREVENTIVES</b>			
PERSONAL PROTECTIVE EQUIPMENT/MATÉRIEL PERSONNEL DE PROTECTION			
GLOVES (SPECIFY)/GANTS (PRÉCISER): <b>Users with special skin problems, cuts or abrasions should wear protective gloves (Latex or PVC).</b>	RESPIRATOR (SPECIFY)/APPAREIL RESPIRATOIRE (PRÉCISER): <b>When working with large amount wear a mask.</b>	EYE (SPECIFY)/APPAREIL OCULAIRE (PRÉCISER): <b>Protective Goggles.</b>	
FOOTWEAR (SPECIFY)/CHAUSSURES (PRÉCISER) <b>N.A.</b>	CLOTHING (SPECIFY)/VÊTEMENTS (PRÉCISER): <b>Use protective clothing.</b>	OTHER (SPECIFY)/AUTRES (PRÉCISER): <b>Do not inhale vapours.</b>	
ENGINEERING CONTROLS (SPECIFY, EG, VENTILATION, ENCLOSED PROCESS)/INSTALLATION TECHNIQUES <b>Use ventilation, inhalation of amine vapours may be irritant to the respiratory system.</b>			
LEAK AND SPILL PROCEDURE/PROCÉDURES EN CAS DE FUITES ET DE DEVERSEMENTS <b>N.A.</b>			
WASTE DISPOSAL/ÉLIMINATION DES RÉSIDUS: <b>Do not empty into drains, water or ground. Dispose of according to local regulations.</b>			
HANDLING PROCEDURES AND EQUIPMENT/MÉTHODES ET ÉQUIPEMENT POUR LA MANUTENTION: <b>Refrigerated storage prolongs shelf life.</b>			
STORAGE REQUIREMENTS/EXIGENCES EN MATIÈRE D'ENTREPOSAGE: <b>Do not store near high temperature or ignition sources. Store away from food and beverages. Avoid high humidity.</b>			
SPECIAL SHIPPING INFORMATION/RENSEIGNEMENTS SPÉCIAUX EN MATIÈRE D'EXPÉDITION <b>Not Regulated.</b>			
<b>SECTION 8 - FIRST AID MEASURES/PREMIER SOINS</b>			
SPECIFIC MEASURES/PREMIERS SOINS PARTICULIERS À ADMINISTRER <b>Eyes: Flush immediately with water for minimum of 15 minutes and seek medical attention.</b> <b>Skin: (Irritation). Wash with soap and water, remove contaminated clothing.</b> <b>Inhalation: Get to fresh air.</b>			
<b>SECTION 9 - PREPARATION DATE OF MSDS/RENSEIGNEMENTS SUR LA PRÉPARATION DE LA FS</b>			
PREPARED BY (GROUP, DEPARTMENT, ETC.)/ PRÉPARÉE PAR (GROUPE, SERVICE, PERSONNE, ETC.)	PHONE NUMBER/NO. DE TÉLÉPHONE	DATE	
 <b>Cadentrix</b>	<b>(416) 994-9356</b>	<b>Rev. 2      1 June 2005</b>	

The statements contained herein are offered for informational purposes only and are based upon technical data that Dentsply Ltd. believes to be accurate. It is intended for use only by persons having the necessary technical skill and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, and FITNESS OR OTHERWISE.