


## DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

SECTION 1 - PRODUCT IDENTIFICATION AND USE/IDENTIFICATION ET UTILISATION DU PRODUIT				
PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUCT		AH 26 Resin		
PRODUCT USE USAGE DU PRODUIT		Dental		
MANUFACTURER'S NAME/NOM DU FABRICANT DENTSPLY DeTray GmbH		SUPPLIER'S NAME/NOM DU FOURNISSEUR DENTSPLY Canada Ltd.		
STREET ADDRESS/NUMERO CIVIQUE De-Tray-Str.1		STREET ADDRESS/NUMERO CIVIQUE 161 Vinyl Court		
CITY/VILLE Konstanz	PROVINCE Germany	CITY/VILLE Woodbridge	PROVINCE Ontario	
POSTAL CODE POSTAL D-78467	EMERGENCY TELEPHONE NO. D'URGENCE 1-49-7531-583-0	POSTAL CODE POSTAL L4L 4A3	EMERGENCY TELEPHONE NO. D'URGENCE (905) 851-6060 or 1-800-263-1437 FAX# (905) 851-9809	
SECTION 2 - HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX				
HAZARDOUS INGREDIENTS/INGREDIENTS DANGEREUX	%	CAS NUMBER/ NUMERO CAS		
Epoxy Resin Product	100	25068-38-6	LD <sub>50</sub> : >5000 mg/Kg	
SECTION 3 - PHYSICAL DATA/CHARACTERISTIQUES PHYSIQUES				
PHYSICAL STATE/ETAT PHYSIQUE Liquid.	ODOUR AND APPEARANCE/ODEUR ET APPARENCE: Yellow liquid.		Chemical Family: Epoxy Resin.	
VAPOUR PRESSURE (mm Hg)/ TENSION DE VAPEUR (mm Hg) <0.01 kPa	VAPOUR DENSITY (AIR-1)/ DENSITE DE LA VAPEUR (AIR-1) N.D.	EVAPORATION RATE/ TAUX D'EVAPORATION N.A.	BOILING POINT (°C)/ POINT D'EBULLITION (°C) N.A.	FREEZING POINT (°C) POINT DE CONGELATION (°C) N.A.
pH: N.A.	SPECIFIC GRAVITY/DENSITE: 1.15-1.20	SOLUBILITY IN WATER: Insoluble.		
SECTION 4 - FIRE AND EXPLOSION DATA/RISQUES D'INCENDIE OU D'EXPLOSION				
FLAMMABILITY/INFLAMMABILITE YES/OUI <input type="checkbox"/> NO/NON <input checked="" type="checkbox"/> IF YES, UNDER WHICH CONDITIONS? SI OUI, DANS QUELLES CONDITIONS?				
MEANS OF EXTINCTION/MOYENS D'EXTINCTION: Foam, CO <sub>2</sub> , Dry chemicals.				
FLASHPOINT (°C) AND METHOD/ POINT D'ECLAIR (°C) ET METHODE DE DETERMINATION 200°C	UPPER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MAXIMAL (% PAR VOLUME) N.A.		LOWER FLAMMABLE LIMIT (% BY VOLUME) SEUIL D'INFLAMMABILITE MINIMAL (% PAR VOLUME) N.A.	
AUTO-IGNITION TEMPERATURE (°C)/ TEMPERATURE D'AUTO-INFLAMMATION (°C): >450°C	HAZARDOUS COMBUSTION PRODUCTS/ PRODUITS DANGEREUX DE COMBUSTION: Poisonous gases may be released during combustion			
SENSITIVITY TO IMPACT/SENSIBILITE A L'IMPACT MECANIQUE EXPLOSION DATA/DONNES SUR LES EXPLOSIONS N.A.		SENSITIVITY TO STATIC DISCHARGE SENSIBILITE A UNE DECHARGE STATIQUE N.A.		
SECTION 5 - REACTIVITY DATA/DONNES SUR LA REACTIVITE				
CHEMICAL STABILITY/STABILITE CHIMIQUE YES/OUI: <input checked="" type="checkbox"/> NO/NON <input type="checkbox"/> IF NO, UNDER WHICH CONDITIONS? SI NON, DANS QUELLES CONDITIONS?				
INCOMPATIBILITY WITH OTHER SUBSTANCES/ INCOMPATIBILITE AVEC D'AUTRES PRODUITS YES/OUI: <input type="checkbox"/> NO/NON <input checked="" type="checkbox"/> IF SO, WHICH ONES? SI OUI, LESQUELS?				
REACTIVITY, AND UNDER WHAT CONDITIONS?/CONDITIONS DE REACTIVITE				
HAZARDOUS DECOMPOSITION PRODUCTS/PRODUITS DE DECOMPOSITION DANGEREUX CO and CO <sub>2</sub>				

## DENTSPLY CANADA LTD. MATERIAL SAFETY DATA SHEET/FICHE SIGNALÉTIQUE

PRODUCT IDENTIFIER/ NOM/IDENTIFICATEUR DU PRODUIT				<b>AH 26 Resin</b>			
<b>SECTION 6 - TOXICOLOGICAL PROPERTIES/PROPRIÉTÉS TOXICOLOGIQUES</b>							
ROUTE OF ENTRY/VOIES D'ABSORPTION							
SKIN CONTACT/ CONTACT DERMIQUE		SKIN ABSORPTION/ ABSORPTION DERMIQUE		EYE CONTACT/ CONTACT OCULAIRE		INHALATION	INGESTION
<input checked="" type="checkbox"/>		<input type="checkbox"/>		<input checked="" type="checkbox"/>		<input type="checkbox"/>	<input checked="" type="checkbox"/>
EFFECTS OF ACUTE EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION AIGUE AU PRODUIT							
N.A.							
EFFECTS OF CHRONIC EXPOSURE TO PRODUCT/EFFETS DE L'EXPOSITION CHRONIQUE AU PRODUIT							
N.A.							
EXPOSURE LIMITS/ SEUILS LIMITES D'EXPOSITION		IRRITANCY OF PRODUCT/ PROPRIÉTÉ IRRITANTE		SENSITIZATION TO PRODUCT/ SENSIBILISATION AU PRODUIT		CARCINOGENICITY/ CANCÉROGENICITÉ	
N.A.		Yes		Yes		N.A.	
TERATOGENICITY/ TERATOGENICITÉ		REPRODUCTIVE TOXICITY/ EFFETS SUR LA REPRODUCTION		MUTAGENICITY/ MUTAGÉNICITÉ		SYNERGISTIC PRODUCTS/ PRODUITS SYNERGIQUES	
N.A.		N.A.		N.A.		N.A.	
<b>SECTION 7 - PREVENTIVE MEASURES/MESURES PREVENTIVES</b>							
PERSONAL PROTECTIVE EQUIPMENT/MATÉRIEL PERSONNEL DE PROTECTION							
<b>Use protective gloves and safety glasses.</b>							
GLOVES (SPECIFY)/GANTS (PRÉCISER):		RESPIRATOR (SPECIFY)/APPAREIL RESPIRATOIRE (PRÉCISER)		EYE (SPECIFY)/APPAREIL OCULAIRE (PRÉCISER):			
<b>Impervious Gloves.</b>		N.A.		<b>Protective Goggles.</b>			
FOOTWEAR (SPECIFY)/CHAUSSURES (PRÉCISER)		CLOTHING (SPECIFY)/VÊTEMENTS (PRÉCISER)		OTHER (SPECIFY)/AUTRES (PRÉCISER)			
N.A.		<b>Protective clothing.</b>					
ENGINEERING CONTROLS (SPECIFY, EG, VENTILATION, ENCLOSED PROCESS)/INSTALLATION TECHNIQUES							
<b>Ensure good ventilation and local exhaust.</b>							
LEAK AND SPILL PROCEDURE/PROCÉDURES EN CAS DE FUITES ET DE DEVERSEMENTS							
<b>Eliminate ignition sources.</b>							
WASTE DISPOSAL/ÉLIMINATION DES RÉSIDUS:							
<b>Disposal must be made according to the local regulations. Do not dispose in municipal waste. Do not allow to reach sewage system. Absorb with liquid-binding material (sand, diatomite, acid binders, sawdust.</b>							
HANDLING PROCEDURES AND EQUIPMENT/MÉTHODES ET ÉQUIPEMENT POUR LA MANUTENTION:							
<b>Store at ambient temperatures out of direct sunlight.</b>							
STORAGE REQUIREMENTS/EXIGENCES EN MATIÈRE D'ENTREPOSAGE:							
<b>Do not allow to reach sewer system or water course. Ensure good ventilation.</b>							
SPECIAL SHIPPING INFORMATION/RENSEIGNEMENTS SPÉCIAUX EN MATIÈRE D'EXPÉDITION							
Not Regulated.							
<b>SECTION 8 - FIRST AID MEASURES/PREMIER SOINS</b>							
SPECIFIC MEASURES/PREMIERS SOINS PARTICULIERS À ADMINISTRER							
<b>Eyes:</b> Flush immediately with water for minimum of 15 minutes and seek medical attention.							
<b>Skin:</b> (Irritation). Wash with soap and water, remove contaminated clothing.							
<b>Ingestion:</b> Induce vomiting and consult physician immediately.							
<b>SECTION 9 - PREPARATION DATE OF MSDS/RENSEIGNEMENTS SUR LA PRÉPARATION DE LA FS</b>							
PREPARED BY (GROUP, DEPARTMENT, ETC.) PRÉPARÉE PAR (GROUPE, SERVICE, PERSONNE, ETC.)		PHONE NUMBER/NO. DE TÉLÉPHONE			DATE		
 Cadentrix		(416) 994-9356			Rev. 2      1 June 2005		

The statements contained herein are offered for informational purposes only and are based upon technical data that Dentsply Ltd. believes to be accurate. It is intended for use only by persons having the necessary technical skill and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, and FITNESS OR OTHERWISE.